



SAFETY CABINETS FOR DANGEROUS PRODUCTS



Cabinets available
in 4 different
colors and 4 different
and in 4 sizes

COMPLIANCE

- Complies with the recommendations of the I.N.R.S.
- Standardized warning labels according to ISO 3864 and ISO 7010 standards, and with European Directive 92/58/EEC.

ADVANTAGES

- Allows the storage of your hazardous products in complete safety and in compliance with legislation.
- Economical solution for storing your hazardous products.
- Supplied with 4 warning labels: «Corrosive» (x2), «Toxic» and «Flammable» (to be affixed according to your use).
- Work surface of the counter cabinets (ref. AZ50 and AZ110) covered by a rubber mat.

ACTIVE SAFETY

- Top and bottom ventilation holes in the doors.
- Hinged door(s) with 110° opening angle.
- 3 points key lock.

PASSIVE SAFETY

- One-piece steel construction.
- Epoxy paint:
- Body grey RAL 7035; Door(s) yellow RAL 1023.
- Body grey RAL 7035; Door(s) red RAL 3000.
- Body grey RAL 7035; Door(s) blue RAL 5010.
- Body grey RAL 7035; Door(s) green RAL 6011.
- Shelves / containment sumps 50 mm high, height-adjustable on racks in 15 mm increments.
- Load capacity per level, evenly distributed, max. 80 kg.
- maximum per level.
- Standardized warning labels for «Flammable», «Corrosive» and «Toxic» products.

OPTIONS

- Additional shelf (ref.: ESE150 and ESE300).
- Galvanized perforated cover plate for shelves / containment sumps (ref.: CE150 and CE300).
- Drilling for connection to a forced ventilation system (ref.: AZVENT), supplied with a Ø 100 mm sleeve (must be ordered with the cabinet).

P/N				Description	External dim. H x W x D (mm)	Storage capacity (L)	Containment capacity (L)	Number of shelves	Weight (kg)
Yellow	Red	Blue	Green						
AZ50	AZ50R	AZ50B	AZ50V	Counter cabinet with 1 door	Ext. : 1000 x 500 x 500 Int. : 890 x 498 x 453	80	20	2	26
AZ110	AZ110R	AZ110B	AZ110V	Counter cabinet with 2 doors	Ext. : 1000 x 930 x 500 Int. : 890 x 925 x 453	155	40	2	43
AZ150	AZ150R	AZ150B	AZ150V	Tall cabinet 1 door	Ext. : 1950 x 500 x 500 Int. : 1840 x 498 x 453	150	40	4	55
AZ300	AZ300R	AZ300B	AZ300V	Tall cabinet 2 doors	Ext. : 1950 x 930 x 500 Int. : 1840 x 925 x 453	300	80	4	73
ESE150	Shelf / base through for AZ50 and AZ150				50 x 488 x 443	-	10	-	3
CE150	Perforated cover plate for AZ50 and AZ150				50 x 483 x 438	-	-	-	5
ESE300	Shelf / base through for AZ110 and AZ300				50 x 915 x 443	-	20	-	6
CE300	Perforated cover plate for AZ110 and AZ300				50 x 913 x 438	-	-	-	10
AZVENT	Hole for ventilation connection				-	-	-	-	-

Refer to the table of chemical incompatibilities and precautions for use, to validate the configuration of your cabinet according to your needs.