



# SAFETY CABINETS FOR DANGEROUS PRODUCTS



◀ AL55



◀ AL115



▲ AL155



▲ AL305J

## COMPLIANCE

- In accordance to articles R.5132-66 and R.5132-68 of the Public Health Code.
- Complies with the recommendations of the I.N.R.S.
- Standardized warning labels in accordance to NF X08.003, ISO 3864 and ISO 7010 standards, and to European Directive 92/58/EEC.

## ADVANTAGES

- Allows the storage of your hazardous products in complete safety and in compliance with the regulations in force.
- Economical solution for the storage of your hazardous products.
- Supplied with 4 Standardized warning labels: «flammable» (x1), «toxic» (x1) and «corrosive» (x2).

## PASSIVE SAFETY

- 8/10th steel construction.
- Electrostatic powder coating in grey RAL 7035 (or yellow RAL 1023 only for the AL300J model).
- Steel retention shelf, 50mm high, adjustable in height on a rack in 25mm increments.

- Steel containment sump, 50mm high at the bottom of the cabinet.
- Standardized warning labels: «flammable» (x1), «toxic» (x1) and «corrosive» (x2).

## ACTIVE SAFETY

- Top and bottom ventilation holes on the sides.
- Hinged door(s) with key lock.

## OPTIONS

- Additional steel shelf / containment sump.
- PVC over-shelf / containment sump.
- Absorbent mat.

**We remind you that these cabinets are not suitable for storing highly corrosive products (nitric, sulphuric, hydrochloric acids, etc.). Please contact us for the range of cabinets best suited to your needs.**

P/N	Description	External dim. H x W x D (mm)	Internal dim. H x W x D (mm)	Storage capacity (L)	Containment capacity (L)	Number of shelves	Load capacity per shelf (kg)	Weight (kg)
AL55	1 door counter cabinet	1000 x 500 x 420	930 x 480 x 390	70	20	2	30	26
AL115	2 door counter cabinet	1000 x 920 x 420	930 x 900 x 390	135	40	2	60	40
AL155	Tall cabinet 1 door	1800 x 500 x 420	1730 x 480 x 390	130	40	4	30	41
AL305	Tall cabinet 2 doors	1800 x 920 x 420	1730 x 900 x 390	240	80	4	60	60
AL305J	Tall yellow cabinet 2 doors	1800 x 920 x 420	1730 x 900 x 390	240	80	4	60	60
ESAL155	Containment shelf / sump for 1 door cabinet	65 x 480 x 375	-	-	10	-	-	3
ESAL305	Containment shelf / sump for 2 doors cabinet	65 x 900 x 375	-	-	20	-	-	6
SEAL155	Containment shelf / sump in PVC for 1 door cabinet	-	-	-	10	-	-	2
SEAL305	Containment shelf / sump in PVC for 2 doors cabinet	-	-	-	20	-	-	3
CAL155	Perforated cover plate for shelf / sump for 1 door cabinet	-	-	-	-	-	-	3
CAL305	Perforated cover plate for shelf / sump for 2 doors cabinet	-	-	-	-	-	-	6
JF10D	Absorbent mats, small and large widths (set of 10 sheets)	-	-	-	-	-	-	1

Refer to the table of chemical incompatibilities and precautions for use, to validate the configuration of your cabinet according to your needs.

Fume hoods, filtering cabinets - Ventilation  
Lithium-ion cabinets  
Flammable cabinets  
Corrosive cabinets  
Toxic cabinets and Multirisk  
Containment and cans  
File cabinets  
Anti-fire equipments  
Showers and first-aid equipments