

# ACTIVE Charcoal FILTERS SELECTION GUIDE

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk














Pesticide range

Fume hoods, filtering cabinets - Ventilation

Containment and Cans

Anti-fire equipments and File cabinets









Showers and first-aid equipments

	Kind of filter	Filter for vapors	P/N from 2011	Weight (kg)	For	P/N before 2011
<b>Filtration boxes and Filtering ventilation cabinets</b>						
  H 150 x Ø 150 (mm)	Organic		CDFORG	3	Range V - Boxes : CDF-A, CDF-P, CVAT, CMAT, CTAT	CDF-F20
	Corrosive		CDFCORG	3		CDF-F35
	Polyvalent <sup>(1)</sup>					CDF-F15
	Formaldehyde		CDDFOR	3		-
	Ammonia		CDFAMM	3		-
  H 70 x W 390 x D 390 (mm)	Organic		ORG50	8	Range H - Boxes : H50	ORG50
	Polyvalent <sup>(1)</sup>		CORG51	9		CORG51
	Formaldehyde		FOR50	8		-
	Ammonia		AMM50	8		-
  H 150 x Ø 150 (mm)	Organic		CDFORG	4	Range 12 Cabinets: AF10	CDF-F20
	Corrosive		CDFCORG	4		CDF-F35
	Polyvalent <sup>(1)</sup>					CDF-F15
	Formaldehyde		CDDFOR	4		-
	Ammonia		CDFAMM	4		-
  H 70 x W 390 x D 390 (mm)	Organic		ORG50	8	Range 12 - Cabinets: AF1, AF2 (before February 2002)	CFR720
	Corrosive		CORG51	9		CFR735
	Polyvalent <sup>(1)</sup>					CFR - CFR70 - CFR715
	Formaldehyde		FOR50	8		-
	Ammonia		AMM50	8		-
 H 20 x W 390 x D 390 (mm)	Polyvalent <sup>(1)</sup>		FCR	1	Range 12 - Cabinets: AF1, AF2 (before February 2002)	FCR
	  H 200 x W 100 x D 400 (mm)	Organic		ORGFC	10	Range 12 - Cabinets: AF1, AF2, AF2-2, AF4-2, AF4-3, AF4-4 (before July 2008) H40
Corrosive			CORGFC	10	FC35	
Polyvalent <sup>(1)</sup>					FC15	
Formaldehyde			FORFC	10	-	
Ammonia			AMMFC	10	-	
  H 70 x W 730 x D 380 (mm)	Organic		ORG300	16	Range 12 - Cabinets: AF1, AF2, AF2-2, AF4-2, AF4-3, AF4-4 (before July 2008)	FC20
	Corrosive		CORG301	15		FC35
	Polyvalent <sup>(1)</sup>					FC15
	Formaldehyde		FOR300	16		-
	Ammonia		AMM300	16		-

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 04/2017

<sup>(1)</sup> organic + corrosive. CORG polyvalent active charcoal filters can adsorb the most common chemical products used: organic and inorganic acids, organic compounds or solvents.

# ACTIVE Charcoal FILTERS SELECTION GUIDE

Models	Kind of filter	Filter for vapors	P/N from 2011	Weight (kg)	For	P/N before 2011
<b>Fume Hoods</b>						
	 H 70 x W 390 x D 390 (mm)	Organic	ORG50	8	Range H Fume hoods : H60 and H70	-
		Polyvalent <sup>(1)</sup>	CORG51	9		-
		Formaldehyde	FOR50	8		-
		Ammonia	AMM50	8		-
	 H 70 x W 545 x D 380 (mm)	Organic	ORG200	11	Range H Fume hoods : H200, H201 H210, H211 PMSF200 H06, H09, H12, H15	ORG200
		Corrosive	CORG201	11,5		COR200
		Polyvalent <sup>(1)</sup>	-	-		CORG201
		Formaldehyde	FOR200	12		-
		Ammonia	AMM200	12		-
		Organic	ORG300 ORG400 ORG1300 ORG2300	16		-
 H 70 x W 730 x D 380 (mm)	Corrosive	Polyvalent <sup>(1)</sup>	CORG301 CORG401 CORG1301 CORG2301	15	Range H Fume hoods : H300, H301, H310, H311, H400, H401, H410, H411, H1300, H1301, H1310, H1311, H2300, H2301, H2310, H2311, PMSF400, PMSF1300, PMSF2300	COR300
						COR400
						COR1300
						COR2300
						COR301
						COR401
Formaldehyde	FOR300 FOR400 FOR1300 FOR2300	16	-			
			-			
			-			
Ammonia	AMM300 AMM400 AMM1300 AMM2300	16	-			
			-			
			-			
 (2x) H 70 x W 545 x D 380 (mm)	Organic	ORG500 ORG1500 ORG2500 ORG3500	22	Range H Fume hoods : H500, H501, H510, H511, H1500, H1501, H1510, H1511, H2500, H2501, H2510, H2511, H3500, H3501, H3510, H3511, PMSF1500	ORG500	
					ORG1500	
					ORG2500	
					ORG3500	
	Corrosive	CORG501 CORG1501 CORG2501 CORG3501	23	COR500		
				COR1500		
	Polyvalent <sup>(1)</sup>	-	-	COR2500		
				COR3500		
	Formaldehyde	FOR500 FOR1500 FOR2500 FOR3500	24	CORG501		
				CORG1501		
Ammonia	AMM500 AMM1500 AMM2500 AMM3500	24	CORG2501			
			CORG3501			
- 806 NU - Filtrair 623/633 MF - Filtrair 936/936-1 MF - Filtrair XL 1046/1346/1646 MF - Filtrair XL 9156 MF - Midcap 633 MF	 H 100 x L 565 x P 405 (mm)	Organic	ORG203	13	Compatible MF3-AS	ORG203
		Corrosive	-	-	Compatible MF3-BE	CORG203
		Polyvalent <sup>(1)</sup>	CORG203	13	Compatible MF3-F	-
		Formaldehyde	FOR203	13	Compatible MF3-K	-
		Ammonia	AMM203	13	Compatible MF4-AS	ORG303
		Organic	ORG303	15	Compatible MF4-BE	CORG303
- Filtrair 814/824/834 HEPA/Carbon - Filtrair XL 1044/1344 MF - Filtrair XL 8104/8124 MF - Midcap 804 MF - Store AVP 804 MF - Store AVP 804 Midcap MF - Store AVPD 804 MF - Store AVPD 804 Midcap MF - Store AVPS 804 MF - Store AVPS Midcap 804 MF - Store AVPSD 804 MF - Store AVPSD Midcap 804 MF - Store AVVE 804 MF	 H 100 x L 755 x P 405 (mm)	Corrosive	CORG303	16	Compatible MF4-F	-
		Polyvalent <sup>(1)</sup>	-	-	-	
		Formaldehyde	FOR303	16	Compatible MF4-K	-
		Ammonia	AMM303	16	-	
		Organic	ORG53	7	Compatible F2BE	-
		Organic	ORG53	6	Compatible F2AS	-
Don't forget, we are able to provide specialised filters for specific applications for fume hoods, filtering ventilation cabinets, filtration boxes... in standard manufacturing and in specific manufacturing, contact-us for more information.	H 404 x L 380 x P 50 (mm)	Polyvalent <sup>(1)</sup>	CORG13	4	Compatible F1BE	-
	H 404 x L 191 x P 50 (mm)	Organic	ORG13	3	Compatible F1AS	-

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 08/2019

(1) organic + corrosive. CORG polyvalent active charcoal filters can adsorb the most common chemical products used : organic and inorganic acids, organic compounds or solvents. We can offer filters for specific applications for all models of fume hoods, air recirculating fume cupboards, safety cabinets... in manufacturing stand or special manufacturing.

- Flammable cabinets
- Corrosive cabinets
- Toxic cabinets and multirisk
- Pesticide range
- Fume hoods, filtering cabinets - Ventilation
- Containment and Cans
- Anti-fire equipments and File cabinets
- Showers and first-aid equipments