













SELECTION GUIDE

You store chemical products		We propose you	Range	You are in compliance with	Pages	
flammables 	Safety cabinets with fire resistance	Cabinets type 105 min	7.90 - 7.90T	EN 14470-1 FM 6050	12, 14	
		Cabinets type 60 min	7.60		18	
		Cabinets type 30 min	3030 - 7030	20, 27		
		FM cabinets	2	FM 6050 – NFPA – OSHA	34	
	Safety cabinets	Safety cabinets with extinguisher	1	-	33	
		Safety storage cabinets	3/5	NF EN ISO 1182	36	
		Secure cabinets	14.E	-	47	
	Safety shelvings	Safety shelvings	R	-	94	
	Containment tanks	Containment tanks	BR	-	91, 93	
			G	EPA – UFC	92	
Safety cans	Stainless steel cans	S	FM – UL	98		
	Steel cans	J	FM – NFPA – OSHA	96		
corrosives 	Safety cabinets built in steel	FM cabinets	4	FM – NFPA – OSHA	41	
	High-corrosive products safety cabinets	Corrosion resistant cabinets	8.A & 8.B	-	42, 43	
		Cabinets type 30 min	3030.M	EN 14470-1	44	
	Safety shelvings	Safety shelvings	R	-	94	
	Containment tanks	Containment tanks	BR	-	91, 93	
			G	EPA – UFC	92	
Safety cans	HPDE cans	J	FM	97		
toxics / noxious 	Safety cabinets	Ventilation cabinets	13	-	53	
		Filtering ventilation safety cabinets	12.X	NF X 15-211*	78	
		Filtering ventilation secure cabinets	14.X		80	
		Filtering ventilation safety cabinets type 30 minutes	3030.X	NF X 15-211* - EN 14470-1	82	
	Safety shelvings	Safety shelvings	R	-	94	
			BR	-	91, 93	
Containment tanks	Containment tanks	G	EPA – UFC	92		
different chemical families (flammables, corrosives, toxics, noxious...) 	Safety cabinets	Under-bench cabinets	6	-	56	
		Multirisk and moduable secure cabinets	14		48	
		Secure cabinets	14.L		50	
		Secure cabinets with drawers	14.T		51	
		Multirisk cabinets type 105 min	Secure cabinets	14.E	47	
			Multirisk cabinets type 90, 60 min	7.90M	16	
			Multirisk cabinets type 30 min	11.B	EN 14470-1	54
			Filtering ventilation safety cabinets	3030.M	24	
			Filtering ventilation safety cabinets	12.X	NF X 15-211*	78
			Filtering ventilation secure cabinets	14.X		80
			Filtering ventilation safety cabinets type 30 minutes	3030.X	NF X 15-211* - EN 14470-1	82
pesticides 	Safety cabinets	Safety cabinets	16.E	-	60	
		Secure cabinets	16.L	-	59	
		Ready to build secure cabinets	16.K	-	60	
	Safety shelvings	Safety shelvings	HPDE safety shelvings	R	EPA – UFC	94
			Steel safety shelvings	R		94
	Containment tanks	Containment tanks	BR	-	91, 93	
G			-	92		
gas cylinders 	Safety cabinets	Antifire safety cabinets	7630BG	EN 14470-2	28	
filing/ document 	Safe for documents	Fire resistant filing/document cabinets	AD73-AD76-AD79	NT Fire 017	108 to 110	
		Safe for documents	AD3 & AD4	NF EN ISO 1182	106,107	
You need to collect your wastes						
flammables 	Safety cabinets	Cabinets type 105 min for recovery of waste	7.90ARD	EN 14470-1	15	
		Cabinets type 30 min for recovery of waste	3030.ARD	EN 14470-1 - FM 6050	26	
You use chemical products						
flammables, corrosives, toxics, noxious... 	Ductless fume hoods	Labopur ductless fume hoods	H	NF X 15-211 – Classe 2	62 to 70, 76	
		PMSF – Ductless fume hood and storage cabinets with recirculating air	H		73	
	Laminar flow clean benches	Laminar flow clean benches	H	EN 12469 - ISO 14644-1 classe 3	74	
	Biological safety cabinets	Biological safety cabinets	H		75	
	Fume hoods with external connection	Fume hoods with external connection	H	-	72	
You would like to vent your cabinet						
ventilation 	Without active charcoal filter	Ventilation boxes	V	-	86	
		ATEX ventilation boxes	V	ATEX	87	
	With active charcoal filter	Recycling air box Labopur	H	NF X 15-211*	84, 85	
		Filtration boxes	V	-	86	
You need equipments						
fire protection 	Safety waste bins	Safety waste bins	P & J	TNO – IBBC – NVBB – FM	103,104	
	Fire blankets	Fire blankets	A	EN 1869-1997	104	
	Safety ashtrays	Safety ashtrays	C	-	105	
first-aid 	Safety combination units, showers and eye-washers	Combination units	DL	EN 15154-1, EN 15154-2	114, 115, 116	
		Showers connected to water network	DL	EN 15154-1	120, 121	
		Eye-washers connected to water network	DL	EN 15154-2	118, 119	
		Portable shower stations	DL	-	123	
		Autonomous eye-washers	DL	-	123	
	Medicine cabinets	Medicine cabinets	19	-	124	

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 04/2017

Flammable cabinets
Corrosive cabinets
Toxic cabinets and multirisk
Pesticide range
Fume hoods, filtering cabinets - Ventilation
Containment and Cans
Anti-fire equipments and File cabinets
Showers and First-aid equipments

*Active coal filters absorption's capacity and filtration system are tested and certified NF X 15-211 - may 2009.