
















ACTIVE CHARCOAL FILTERS SELECTION GUIDE

MODELS	FILTER TYPES	VAPOR FILTERS	P/N FROM 2021	WEIGHT (KG)	FOR	P/N BEFORE 2020
 	 H 70 x L 545 x P 380 (mm)	Organic	ORG200	11	Fume hoods Range H : H06+, H09+, H12+, H15+, H18+, H200, PMSF200	ORG500 ⁽²⁾ ORG1500 ORG2500 ORG3500
		Corrosive	CORG201	11,5		COR500 ⁽²⁾ COR1500 COR2500 COR3500
		Polyvalent ⁽¹⁾			Cabinets Range 12+ : AF2+, AF22+, AF42+, AF43+, AF44+,	CORG501 ⁽²⁾ CORG1501 CORG2501 CORG3501
		Formaldehyde	FOR200	12		Two filters required for P/N : H500, H1500, H2500, H3500, PMSF500, PMSF1500, PMSF2500, PMSF3500
		Ammonia	AMM200	12	AMM500 ⁽²⁾ AMM1500 AMM2500 AMM3500	
	 H 70 x L 730 x P380 (mm)	Organic	ORG300	16	Fume hoods Range H : H300, H400, H1300, H2300, PMSF300, PMSF400, PMSF1300, PMSF2300	ORG400 ORG1300 ORG2300
		Corrosive	CORG301	15		COR300 COR400 COR1300 COR2300
		Polyvalent ⁽¹⁾			Cabinets range 12X : AF2X, AF2-2X, AF4-2X, AF4-3X, AF4-4X, AF1X, AF1X1, AF2X1, AF2C, AF22C, AF42C, AF43C, AF44C, AF1C, AF1C1 et AF2C1. (After July 2008)	CORG401 CORG1301 CORG2301
		Formaldehyde	FOR300	16		FOR400 FOR1300 FOR2300
		Ammonia	AMM300	16	AMM400 AMM1300 AMM 2300	
- 806 NU - Filtair 623/633 MF - Filtair 936/936-1 MF - Filtair XL 1046/1346/1646 MF - Filtair XL 9156 MF - Midcap 633 MF	 H 100 x L 565 x P 405 (mm)	Organic	ORG203	13	Compatible MF3-AS	ORG203
		Corrosive	CORG203	13	Compatible MF3-BE	CORG203
		Polyvalent ⁽¹⁾				
		Formaldehyde	FOR203	13	Compatible MF3-F	
		Ammonia	AMM203	13	Compatible MF3-K	
- Filtair 814/824/834 HEPA/Charcoal - Filtair XL 1044/1344 MF - Filtair XL 8104/8124 MF - Midcap 804 MF - Store AVP 804 MF - Store AVP 804 Midcap MF - Store AVPD 804 MF - Store AVPD 804 Midcap MF - Store AVPS 804 MF - Store AVPS Midcap 804 MF - Store AVPSD 804 MF - Store AVPSD Midcap 804 MF - Store AVVE 804 MF	 H 100 x L 755 x P 405 (mm)	Organic	ORG303	15	Compatible MF4-AS	ORG303
		Corrosive	CORG303	16	Compatible MF4-BE	CORG303
		Polyvalent ⁽¹⁾				
		Formaldehyde	FOR303	16	Compatible MF4-F	-
		Ammonia	AMM303	16	Compatible MF4-K	-
- Schoolcap AVPS 802 - Store 832	H 404 x L 380 x P 50 (mm)	Polyvalent ⁽¹⁾	CORG53	7	Compatible F2BE	-
		Organic	ORG53	6	Compatible F2AS	-
		Polyvalent ⁽¹⁾	CORG13	4	Compatible F1BE	-
- Ministore 822 - Shelf 812	H 404 x L 191 x P 50 (mm)	Organic	ORG13	3	Compatible F1AS	-

© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 10/2023

MODELS	FILTER TYPES	VAPOR FILTERS	P/N FROM 2021	WEIGHT (KG)	FOR	P/N BEFORE 2020
	 H 150 x Ø 150 (mm)	Organic	CDFORG	4	Box Range V: CDF-A, CDF-P, CVAT, CMAT, CTAT Cabinet Range 12 : AF10	CDF-F20
		Corrosive	CDFCORG	4		CDF-F35
		Polyvalent ⁽¹⁾				CDF-F15
		Formaldehyde	CDFFOR	4		-
		Ammonia	CDFAMM	4		-
  	 H 70 x L 390 x P 390 (mm)	Organic	ORG50	8	Fume hoods Range H : H70 Box Range H : H50 Cabinets Range 12 : AF1, AF2 (before February 2002)	CDF-F720
		Corrosive	CORG51	9		CDF735
		Polyvalent ⁽¹⁾				CFR CFR70 CFR715 CORG51
		Formaldehyde	FOR50	8		-
		Ammonia	AMM50	8		-
	 H 200 x L 100 x P 400 (mm)	Organic	ORGFC	10	Cabinets range 12 : AF1, AF2, AF2-2, AF4-2, AF4-3, AF4-4 (before July 2008) H40	FC20
		Corrosive	CORGFC	10		FC35
		Polyvalent ⁽¹⁾				FC15
		Formaldehyde	FORFC	10		-
		Ammonia	AMMFC	10		-

Fume hoods, filtering cabinets - Ventilation

Lithium-ion cabinets

Flammable cabinets

Corrosive cabinets

Toxic cabinets and Multirisk

Containment and cans

File cabinets
Anti-fire equipments

Showers and first-aid equipments

(1) CORG multipurpose active charcoal filters can adsorb the vast majority of commonly used products: corrosive products, organic compounds or solvents.
 We can offer you filters for specific applications for all models of fume hoods, recirculating fume hoods, safety cabinets, etc. in standard or special production.

+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+

© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations.10/2023