

SECURIFLAMM SAFETY WASTE BINS



▲ 50RN

▲ 30MN

▲ 18GP

COMPLIANCE

- Compliant to CNPP FSE 8414-D, TNO and IBBC approved No. B 89.623 ANPI and NVBB No. NT 108.

BENEFITS

- In the event of a fire internal the bin, the funnel-shaped smothering lid prevents hot gases from escaping and sends them back to the fire, which, deprived of oxygen, extinguishes itself.
- Many models available (capacity and colors).

TECHNICAL CHARACTERISTICS

- 8/10th steel construction.
- Anti-cut edges on the upper part and anti-scratch on the lower part.
- Funnel-shaped safety lid making the whole unit self-extinguishing.
- Epoxy paint: pearl grey lid and body; black lid and red body; black lid and metal grey body.
- Capacity: 15, 30, 50, 80 and 110 litres.

P/N	Color	Dimensions H x Ø (mm)	Capacity (L)	Weight (kg)
18GP	Pearl grey body and lid	360 x 255	15	1,7
18RN	Red body - Black cover	360 x 255	15	1,7
18MN	Metal grey body - Black lid	360 x 255	15	1,7
30GP	Pearl grey body and lid	470 x 335	30	3,9
30RN	Red body - Black lid	470 x 335	30	3,9
30MN	Metal grey body - Black lid	470 x 335	30	3,9
50GP	Pearl grey body and lid	625 x 335	50	4,6
50RN	Red body - Black lid	625 x 335	50	4,6
50MN	Metal grey body - Black lid	625 x 335	50	4,6
80GP	Pearl grey body and lid	510 x 450	80	8
80RN	Red body - Black lid	510 x 450	80	8
80MN	Metal grey body - Black lid	510 x 450	80	8
110GP	Pearl grey body and lid	710 x 450	110	9,3
110RN	Red body - Black lid	710 x 450	110	9,3
110MN	Metal grey body - Black lid	710 x 450	110	9,3

© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 09/2021

Fume hoods, filtering cabinets - Ventilation

Lithium-ion cabinets

Flammable cabinets

Corrosive cabinets

Toxic cabinets and Multirisk

Containment and cans

File cabinets
Anti-fire equipments

Showers and first-aid equipments