

## SAFETY CABINETS WITH EXTINGUISHERS FOR FLAMMABLE PRODUCTS



▲ AMS10



▲ AMS20



1

CDVA



2

EX500G



▲ AMS100



▲ AMS200

### COMPLIANCE

- Compliant with CNPP FSE 98-006 specifications.
- Standardized warning labels in accordance to ISO 3864, ISO 7010 standards and to the European directive 92/58/EEC.

### ADVANTAGES

- Large containment sump on the AMS100 and AMS200 models.
- CO<sub>2</sub> fire extinguisher with automatic activation in case of fire.
- Sloping roof to prevent unauthorised storage.

### PASSIVE SAFETY

- 15/10th steel construction.
- Epoxy paint as standard:
- Beige RAL 1013 for models AMS10 and AMS20,
- Yellow RAL 1023 for models AMS100 and AMS200.
- Standardized «Flammable» symbol.
- 3 adjustable perforated shelves on racks (except AMS10 and AMS20, one shelf, not height adjustable).
- Containment sump integrated in the lower part with drain plug.

### ACTIVE SAFETY

- Equipped with a CO<sub>2</sub> extinguisher.
- 120 g for the AMS10 and AMS20 models,
- 500 g for models AMS100 and AMS200.
- Automatic triggering of the extinguisher by thermal bulb between 65° C and 72° C.
- Top and bottom side vents Ø 100 mm with flame arrestor grid (except AMS10 and AMS20).
- Double swing doors (except AMS10 and AMS20, one door).
- Self-closing doors.
- Key lock.

### OPTIONS

- Additional shelf for AMS100 and AMS200.
- Spare CO<sub>2</sub> extinguisher.
- Cabinets also available without CO<sub>2</sub> extinguisher.
- Absorbent mats to be placed at the bottom of the Containment sump (packs of 10 sheets, reference JF10D).
- Ventilation or filtration boxes.
- 1 To improve the longevity and safety of your cabinet, we recommend connecting it to a forced ventilation system. For more information, please see the ventilation section of this catalogue and do not hesitate to contact us.

P/N	External dim. H x W x D (mm)	Internal dim. H x W x D (mm)	Delivered with a CO <sub>2</sub> extinguisher	Storage (L)	Retention(L)	Number of shelves	Load capacity per shelf (kg)	Weight (kg)
AMS10	455 x 500 x 250	260 x 490 x 220	120 g	10	10	1	20	12
AMS20	500 x 700 x 250	260 x 695 x 230	120 g	20	15	1	25	20
AMS100	1500 x 800 x 400	1250 x 780 x 370	500 g	100	112	3	50	70
AMS100SE			non					
AMS200	1800 x 1000 x 400	1530 x 980 x 370	500 g	200	130	3	75	101
AMS200SE			non					
ESAMS1	30 x 780 x 370	Additional shelf for AMS100 and AMS100SE					50	8
ESAMS2	30 x 980 x 370	Additional shelf for AMS200 and AMS200SE					75	10
EX120G	Replacement CO <sub>2</sub> extinguisher for AMS10 and AMS20						-	0,2
EX500G	Replacement CO <sub>2</sub> extinguisher for AMS100 and AMS200						-	0,6