



VENTILATED CABINETS FOR EXTERNAL CONNECTION



AF2+XT



AF22+XT



AF42+XT



AF43+XT



AF44+XT

TECHNICAL DATA

Storage compartment(s) :

- 15/10th steel construction, double-walled doors
- Large, lockable polycarbonate glass door(s).
- White RAL 9010 and blue RAL 5015 polyester paint.
- Height-adjustable perforated shelves.

- Roof ventilation Ø 100mm

«Ventilation/filtration» compartment:

Silent fan: CE 220-240 volts - 50 Hz; flow rate 202 m³/hour; 0.2 A; 70 watts; sound level 49 dB(A) (excluding airflow) with plug and cable.

| P/N | AF2+XT | AF22+XT | AF42+XT | AF43+XT | AF44+XT |
|---|-------------------------|----------------------|----------------------|------------------------------|----------------------|
| External dimensions H x W x D (mm) | 2020 x 800 x 650 | 2020 x 800 x 650 | 2020 x 1600 x 650 | 2020 x 1600 x 650 | 2020 x 1600 x 650 |
| Internal dimensions H x W x D (mm) | 1660 x 780 x 500 | 1660 x (2x380) x 500 | 1660 x (2x780) x 500 | 1660 x (1x780 + 2x380) x 500 | 1660 x (4x380) x 500 |
| Compartments | 1 | 2 | 2 | 3 | 4 |
| Shelves | 4 | 8 halves | 8 | 4+8 halves | 16 halves |
| Retention bin(s) | 1 | 2 | 2 | 3 | 4 |
| Storage capacity (L) | with 4 shelves | 140 | 120 | 280 | 260 |
| | with 5 shelves | 180 | 150 | 360 | 360 |
| Weight (kg) | 140 | 150 | 280 | 290 | 300 |
| Connecting element for use without filter | CAR200 (3 kg) (include) | | | | |
| Ø 100mm connection kit | KRC (include) | | | | |

FILTER FANNED PANEL OR WALL CABINET



| P/N | Description | Dimensions H x W x D (mm) | Number of shelves | Number of containment sump | Storage capacity in 1L or 2.5L bottles | Load per shelf | Weight (kg) |
|------|------------------------------------|---|-------------------|----------------------------|--|----------------|-------------|
| AF11 | Cabinet without ventilation system | External 550 x 810 x 325 mm Internal 535 x 790 x 275 | 1 | 1 | 60 or 10 | 20 | 35 |

© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 09/2021

Fume hoods, filtering cabinets - Ventilation

Lithium-ion cabinets

Flammable cabinets

Corrosive cabinets

Toxic cabinets and Multirisk

Containment and cans

File cabinets
Anti-fire equipments

Showers and first-aid equipments