

RANGE 14.T



SECURE CABINETS WITH SLIDING DRAWERS FOR DANGEROUS PRODUCTS



ALT155 ▲



ALT305 ▲



Locking system drawers



Shelves/retention on coulissants rails



Vents



Handle with integrated lock



Fix on the wall 2 points

CERTIFICATION

- Large warning labels in compliance with ISO 3864, NF X08.003 norm and the European directive 92/58/CEE.

ADVANTAGES

- Sliding containment shelves for an easy access to the stored products.
- Locking system for the others drawers in the retracted position when one drawer is already pull-out.
- Possibility to fix on the wall.
- Supplied with "Toxic", "Flammable" or "Acid (2x)" large warning labels following your needs (to stick on according to your storage).

PASSIVE SECURITY

- Built in steel 8/10th.
- White grey paint RAL 7035.
- 3 steel sliding containment shelves (height: 50 mm; maximum load by shelf: 45 kg).
- 1 steel sliding containment sump (height: 80 mm).

ACTIVE SECURITY

- High and low vents.
- Handling door(s) with key locking and integrated handle.

OPTIONS

- Industrial absorbent mat (ref. JF10D).

P/N	Description	External dim. H x W x D (mm)	Internal dim. H x W x D (mm)	Storage capacity (L)	Containment capacity (L)	Number of shelves/retention	Maximum load by shelf ⁽¹⁾	W. (kg)
ALT155	Tall cabinet with drawers 1 door	1950 x 500 x 500	1718 x 379 x 430	150	40	4	30	72
ALT305	Tall cabinet with drawers 2 doors	1950 x 916 x 500	1718 x 795 x 430	300	85	4	50	102
SEALT155	Containment shelf/sump in PVC for cabinet 1 door						-	2
SEALT305	Containment shelf/sump in PVC for cabinet 2 doors						-	3
JF10D	Industrial absorbent mat, small and large widths (pack of 10 sheets)						-	1

⁽¹⁾ Evenly distributed

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 06/2018

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk

Pesticide range

Fume hoods, filtering cabinets - Ventilation

Containment and Cans

Anti-fire equipments and Files cabinets

Showers and first-aid equipments