

ECOSAFE

TRIONYX

138 avenue Henri Barbusse
92140 CLAMART
France
Ph.: +33 (0)1 46 45 80 00
Fax: +33 (0)1 46 45 85 98

Local representative: ECOSAFE

Chemin des Champs-Courbe
1024 ECUBLENS
Suisse
Ph.: +41 (0) 694 10 04
Fax: +41 (0) 21 694 10 00

www.ecosafesa.com

Jean-Michel THIEULIN
sales@ecosafesa.com

www.trionyx.fr

Ecosafe is the swiss specialist in storage and use of dangerous products (flammable, corrosive, toxic, noxious and pesticides products) in laboratory and industry. We offer a large range of safety cabinets for the storage of chemical products (fire resistant cabinets, corrosive resistant cabinets, ventilated cabinets...) and a range of ductless fume hoods equipped with active coal filters or particles filters (HEPA H14). We also propose a range of safety showers and eye-washers, antifire equipments et pallet and accumulation centers for spill containment. Our materials are in compliance with international regulations and norms: EN 14470-1, EN 14470-2, FM, NF X 15-211. Our products are checked before expedition in compliance with the ISO 9001-2008 norm.

Ecosafe ist Experte für die sichere Lagerung und Handling von Gefahrstoffen (leicht entzündbare, korrosive oder giftige Stoffe und vieles mehr) für Laboratorien und Industrien.

Wir bieten Ihnen eine grosse Auswahl von verschiedenen Sicherheitsschränken für die Lagerung Ihrer chemischen Produkte und Gefahrstoffen an, ebenso wie eine Vielzahl von Umluftfilterabzüge mit Aktivkohlefilter und/ oder Sonderfilter (HEPA H14).

Auch haben wir eine Auswahl an Sicherheitsduschen und Augenspülungen, feuerwiderstandsfähige Ware sowie Produkte, die bestimmte Stoffe zurückhalten bzw. auffangen.

Unsere Ware unterlaufen regelmässige Kontrollen und sind von folgenden Normen zertifiziert: EN 14470-1, EN 14470-2, FM, NF X 15-211.

Alle unsere Produkte werden vor dem Versand unter der Norm ISO 9001-2008 kontrolliert.

