

RANGE 14



MODULABLE MULTIRISK SAFETY CABINETS



AS [▲]

A2PCA [▲]

A3PCA [▲]

A2+2PCA [▲] ▲ ASVM11

▲ ASV



CERTIFICATION

- Large warning labels in compliance with ISO 3864 norm and the European directive 92/58/CEE.

ADVANTAGES

- Possibility to separate the dangerous product families.
- Modulable solution for the storage of your chemicals.
- Inside equipments of the different compartments following your needs.
- Possibility to join the different models together (except AB56 model; only A1T, A2T, A3T and A4T can be joined together).

PASSIVE SECURITY

- Built in steel 10/10th.
- White polyester paint RAL 9010.
- "Toxic", "Flammable" or "Acid" (x2) large warning label following your needs
- Steel containment shelves adjustable in height.
- One removable base trough in steel in the bottom of each compartment.

ACTIVE SECURITY

- High and low vents :
 - Outlet with flame arresters;
 - Ø 100mm outlet for a possible external connection.
- Handling door(s) (except for A1T type, sliding door) with key locking system.

OPTIONS

- PVC equipment for "Acids" compartment including PVC protection shelf and PVC retention tank.
- Additional steel retention shelf.
- Industrial absorbent mat (Ref JF10D).

Ventilated under-bench cabinet with two compartments Model AB56⁽¹⁾

Delivered with one electric fan fixed in the right compartment (230 volts - 50 Hz - 40 Watts - 0,19 A - 54 dB(A), excluding air flow - flow 150 m³/h - protection IP44 - Exhaust outlet Ø 100 mm) and PVC retention shelf and base trough.



AB56 [▲]

N.B.: Do not store flammable products and solvents in the compartment with a fan.

⁽¹⁾Cabinet without fan : ref. AB56SV

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 08/2015

Refer to the table of chemical incompatibilities and precautions for use, to validate the configuration of your cabinet according to your needs.

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk

Pesticide range


Fume hoods, filtering cabinets - Ventilation

Containment and Cans

Anti-fire equipments and Files cabinets

Showers and first-aid equipments

RANGE 14




P/N	ASP	ASDP	AB56	AB56SV
Door(s)	1 Handling	2 Handling	2 Handling	
Compartment(s)	1	1	2	
Dimensions H x W x D (mm)	670 x 560 x 545	670 x 1100 x 545	670 x 1100 x 545	
Storage cap. (L)	30	60	60 (2x 30)	
Containment cap. (L)	30	60	60 (2x 30)	
Containment shelf steel/PVC	1 / option	1 / option	2 / 1	
Steel/PVC containment sump	1 / option	1 / option	1 / 1	
Maximum load by shelf ⁽¹⁾	50	50	50	
Weight (kg)	25	48	50	



^ A1T



P/N	AS	ASV	A2P	A2PCA	A3P	A3PCA	A1T	A2T	A3T	A4T	
Door(s)	1 Handling	1 Glazed	2 Handling		3 Handling		of 1 to 4 Handling with drawers				
Compartment(s)	1	1	2		3		1	2	3	4	
Dimensions H x W x D (mm)	1825 x 560 x 545	1825 x 560 x 545	1825 x 560 x 545		1825 x 560 x 545		1825 x	560 / 1120 / 1680 / 2240			x 545
Storage cap. (L)	150	150	110 (2x 55)		90 (3x 30)		120	240	360	480	
Containment cap. (L)	40	40	70 (2x 35)		90 (3x 30)		20	40	60	80	
Containment shelf steel/PVC	3 / option	3 / option	4 / option	4 / 2	3 / option	3 / 1	4 / option	8 / option	12 / option	16 / option	
Steel/PVC containment sump	1 / option	1 / option	2 / option	1 / 1	3 / option	2 / 1	- / -	- / -	- / -	- / -	
Maximum load by shelf ⁽¹⁾	50	50	50		50		50	50	50	50	
Weight (kg)	55	55	85		95		80	160	240	320	



P/N	ASDD	ASDV	ASDM11	ASVM11	A2+2P	A2+2PCA	ASDM22	ASVM22
Doors	2 Handling	2 Glazed	2 Handling	2 Glazed	4 Handling		2 Handling	2 Glazed
Compartment(s)	1	1	2	2	4		4	4
Dimensions H x W x D (mm)	1825 x 1100 x 545	1825 x 1100 x 545	1825 x 1100 x 545	1825 x 1100 x 545	1825 x 1100 x 545	1825 x 1100 x 545	1825 x 1100 x 545	1825 x 1100 x 545
Storage cap. (L)	300	300	260 (2x 130)	260 (2x 130)	220 (4x 55)		220 (4x 55)	220 (4x 55)
Containment cap. (L)	80	80	92 (2x 46)	92 (2x 46)	140 (4 x 35)		140 (4 x 35)	140 (4 x 35)
Containment shelf steel/PVC	3 / option	3 / option	6 / option	6 / option	8 / option	8 / 2	8 / 2	8 / 2
Steel/PVC contain- ment sump	1 / option	1 / option	2 / options	2 / options	4 / option	3 / 1	3 / 1	3 / 1
Maximum load by shelf ⁽¹⁾	50	50	50	50	50		50	50
Weight (kg)	95	95	115	115	170		130	130

P/N	EE1	EE2	ES1T	ESASDD	SEPVC	SBPVC	EX120G	EX500G
Description	Additional shelf	Additional shelf	Additional shelf	Additional shelf	Protection shelf in PVC	PVC Containment sump	CO ₂ extinguisher	
For cabinet	ASP, AS, ASV, A2P, A3P, A2+2P	ASDM11, ASVM11, ASDM22, ASVM22	A1T, A2T, A3T, A4T	ASDP, ASDD, ASDV	All models except A1T (ref.: ESP1T), for models wide in 1100 to 1 compartment provide 2 pieces per level		Contact us	
Dimensions H x W x D (mm)	22 x 550 x 465	30 x 545 x 468	30 x 420 x 450	20 x 1090 x 465	20 x 530 x 460	90 x 530 x 460	-	
Containment cap (L)	5	7	5	10	5	21	-	
Weight (kg)	3	3	3	6	2	3	0,2	0,6

⁽¹⁾ Evenly distributed

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 09/2019

Flammable
cabinets

Corrosive
cabinets

Toxic cabinets
and multirisk

Pesticide
range

Fume hoods, filtering
cabinets - Ventilation

Containment
and Cans

Anti-fire equipments
and Files cabinets

Showers and first-
aid equipments