





SAFETY CABINETS WITH SLIDING DRAWERS FOR DANGEROUS PRODUCTS




1
Wall mounting
in 2 points




2
Drawer blocking
system



3
Shelves/tubs on
sliding rails



4
Ventilation



5
Handle with
integrated lock

COMPLIANCE

- In accordance with articles R.5132-66 and R.5132-68 of the Public Health Code.
- Complies with the recommendations of the I.N.R.S.
- Standardized warning labels according to NF X08.003, ISO 3864 and ISO 7010 standards, and with European Directive 92/58/EEC.

ADVANTAGES

- Sliding drawers for an easy access to the stored products.
- Drawers are locked in the retracted position when a drawer is already pulled out.
- Possibility of wall mounting.
- Supplied with 4 Standardized warning labels: «flammable» (x1), «toxic» (x1) and «corrosive» (x2).

PASSIVE SAFETY

- 8/10th steel construction.
- Grey RAL 7035 electrostatic powder coating.
- 3 steel drawers, 50mm high, mounted on sliding rails (load capacity per shelf evenly distributed: 30 to 50 kg).
- Steel containment sump at the bottom of the cabinet, 80 mm high, mounted on sliding rails.
- Standardized warning labels: «flammable» (x1), «toxic» (x1) and «corrosive» (x2) to be affixed as required.

ACTIVE SAFETY

- Top and bottom ventilation holes.
- Hinged door(s) with key lock, integrated handle.

OPTIONS

- Additional steel shelf/Containment sump.
- Absorbent matting.

P/N	Description	Dimensions H x W x D (mm)	Internal dim. H x W x D (mm)	Storage capa- city (L)	Containment capacity (L)	Number of shelves	Load capacity per shelf (kg)	Weight (kg)
ALT155	1 door tall cabinet with drawers	1950 x 500 x 500	1718 x 379 x 430	150	40	4	30	72
ALT305	2 doors tall cabinet with drawers	1950 x 916 x 500	1718 x 795 x 430	300	85	4	50	102
SEALT155	Containment shelf / sump in PVC for 1 door cabinet						-	2
SEALT305	Containment shelf / sump in PVC for 2 doors cabinet						-	3
JF10D	Absorbent mats, small and large widths (set of 10 sheets)						-	1

Refer to the table of chemical incompatibilities and precautions for use, to validate the configuration of your cabinet according to your needs.

© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 09/2021

Fume hoods, filtering
cabinets - Ventilation

Lithium-ion
cabinets

Flammable
cabinets

Corrosive
cabinets

Toxic cabinets
and Multirisk

Containment
and cans

File cabinets
Anti-fire equipments

Showers and first-
aid equipments