



Fume hoods, filtering cabinets - Ventilation  
Lithium-ion cabinets  
Flammable cabinets  
Corrosive cabinets  
Toxic cabinets and Multirisk  
Containment and cans  
File cabinets  
Anti-fire equipments  
Showers and first-aid equipments

## SAFETY CABINETS - EN 14470-1 60 MINUTES FOR FLAMMABLE PRODUCTS



▲ 765+E + C235



1 Standardized labeling



2 Self-closing door with thermofuse



3 Key locking system



4 Retention shelves



5 Moving by pallet truck (empty)



6 Removable containment sump

### COMPLIANCE

- Tested and approved to EN 14470-1 (04/2004) and EN 1363-1 (06/2000).
- Standardized warning labels in accordance to ISO 3864, Directive 92/58/EEC and EN 14470-1 (04-2004).

### ADVANTAGES

- 60 minutes fire resistance according to the ISO 834 standard curve.
- Fire tested by a French test laboratory accredited to EN:17025 by COFRAC, according to the ISO 834 standard curve.
- Full height handle for better ergonomics and ease of opening
- Legs for easy movement by pallet truck: avoids heavy handling and allows the cabinet to be moved easily in the event of a new layout
- Aesthetic and easy to use skirting board cover without the need for fixing tools
- Castors included for models 761+E and 762+E.

### PASSIVE SAFETY

- Double wall construction:
- Outer walls in 12/10th grade steel, painted in white epoxy RAL 9010.
- Internal walls in white RAL 9010 melamine.
- Thermal insulation panels between the walls to limit thermal bridges.
- Removable containment sump in the lower part. The containment sump must be positioned at the bottom of your cabinet. If you wish to use it as storage space, you must first place a perforated cover plate (see options) corresponding to your containment sump.

- Standardized labeling
- Adjustable retention shelves (at 6.4 cm intervals).
- Fixing point for earthing.

### ACTIVE SAFETY

- Ventilation holes with Ø 100 mm ventilation outlet for possible connection.
- Ventilation ducts with hot-melt system to isolate the contents of the enclosure in the event of fire.
- Thermo-expanding door seals.
- Self-closing door(s) with key lock.
- Door(s) holding system in open position with automatic closing thermo-regulated at 47°C.

### OPTIONS

- Additional shelf
- Perforated cover plate for containment sump.
- Sliding drawer/Containment sump
- PVC overlay
- Ventilation or filter boxes
- Absorbent mats to be placed at the bottom of the containment sump (packs of 10 sheets)
- Yellow doors available on 764+ and 765+ models
- Code lock (1 per door)

# RANGE 760+



## PRE-EQUIPPED CABINETS

P/N pour porte		Type	Number of door(s)	External dimensions H x W x D (mm)	Internal dimensions H x W x D (mm)	Storage capacity (L)	Delivered with		Maximum load by shelf (kg)	Weight (kg)
white	yellow						Retention shelf	Bac de rétention		
761+E	-	under-bench	1	630 x 635 x 570 <sup>(1)</sup>	420 x 490 x 285	30	1	1	30	127
762+E	-	under-bench	2	630 x 1137 x 570 <sup>(1)</sup>	420 x 1000 x 285	60	1	1	50	213
768+E	-	counter	1	1100 x 635 x 670 <sup>(2)</sup>	820 x 490 x 410	65	2	1	30	182
763+E	-	counter	2	1100 x 1137 x 670 <sup>(2)</sup>	820 x 1000 x 410	115	2	1	50	297
764+E	794+PJE	tall	1	1950 x 635 x 620 <sup>(3)</sup>	1620 x 490 x 410	130	3	1	30	287
765+E	795+PJE	tall	2	1950 x 1137 x 620 <sup>(3)</sup>	1620 x 1000 x 410	250	3	1	50	453
765+900E	-	tall	2	1950 x 905 x 620 <sup>(3)</sup>	1620 x 765 x 410	200	3	1	46	438

## ACCESSORIES

P/N	Description	External dimensions H x W x D (mm)	Weight (kg)
E1	Retention shelf for 1 door cabinet type 761+	30 x 485 x 285	3
E2	Retention shelf for 2 door cabinet type 762+ <sup>(1)</sup>	30 x 992 x 285	3
E48	Retention shelf for 1 door cabinet - W: 635 - D: 670 mm <sup>(1)</sup>	30 x 490 x 410	3
E35	Retention shelf for 2 door cabinet - W: 1135 - D: 670 mm	30 x 992 x 410	4
E235-900	Retention shelf for 2 door cabinet - W: 900 - D: 670 mm	30 x 760 x 410	4
B1	Containment sump for 1 door cabinet type 761+	100 x 480 x 315	3
B2	Containment sump for 1 door cabinet type 762+ - W: 900 - D: 670 mm	100 x 985 x 315	6
B148	Containment sump for 1 door cabinet - W: 635 - D: 670 mm	100 x 480 x 410	4
B235	Containment sump for 2 door cabinet - W: 1135 - D: 670 mm	100 x 985 x 410	8
B235-900	Containment sump for 2 door cabinet - W: 900 - D: 670 mm	100 x 753 x 410	5
C1	Perforated cover plate for 1 door cabinet type 761+	90 x 470 x 305	3
C2	Perforated cover plate for 2 doors cabinet type 762+ -W: 1135 - D: 670 mm	90 x 980 x 305	6
C148	Perforated cover plate for 1 door cabinet - W: 635 - D: 670 mm	90 x 460 x 395	3
C235	Perforated cover plate for 2 doors cabinet - W: 1135 - D: 670 mm	90 x 975 x 395	6
C235-900	Perforated cover plate for 2 doors cabinet - W: 900 - D: 615 mm	90 x 747 x 400	5
SERCODE	Digital code lock	-	1

(1) Depth with connecting sleeve at the rear. Depth without sleeve: 520 mm

(2) Depth with connecting sleeve at the rear. Depth without sleeve: 620 mm

(3) Height with connecting sleeve above. Height without sleeve: 1900 mm

+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+

© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 09/2021

Fume hoods, filtering cabinets - Ventilation  
Lithium-ion cabinets  
Flammable cabinets  
Corrosive cabinets  
Toxic cabinets and Multirisk  
Containment and cans  
File cabinets  
Anti-fire equipments  
Showers and first-aid equipments