




NON-MAGNETIC SAFETY CABINET TALL, 2 GLASS HINGED DOORS



▲ AA220V

CONCEPTION

Specially designed for the storage of equipment sensitive to external magnetic fields

DESCRIPTION

- Cabinet construction in cellular PVC and austenitic stainless steel.
- Colour of doors: white RAL 9010.
- Number of doors: 2.
- 5 shelves, 4 of which are removable and 1 fixed.
- Hinged doors on plastic hinges with austenitic stainless steel axle, manual closing.
- Glass or PMMA front (3 per door, not partially UV filtering).
- Type of lock: 1 austenitic stainless steel 1-point swing lock
- Plated fields: ABS 8/10 mm white.

ADVANTAGES

- Made in Europe.
- Adjustable plastic feet for levelling.
- Fixed and removable shelf in cellular PVC on brackets with reinforcements.
- Construction in cellular PVC and stainless steel: magnetically neutral.
- Honeycomb ventilation holes.

OPTION

- Wall fixing kit.

Réf.	Designation	Ext.dim. H x W x D (mm)	Removable shelves	Fixed shelf	Load capacity per shelf (kg)	Weight (kg)
AA220V	Non-magnetic cabinet with glass front (does not filter UV light)	2015 x 1200 x 800	4	1	40	150
AA220P	Non-magnetic cabinet PMMA front (partially UV filtering)	2015 x 1200 x 800	4	1	40	151

© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 09/2021

- Fume hoods, filtering cabinets - Ventilation
- Lithium-ion cabinets
- Flammable cabinets
- Corrosive cabinets
- Toxic cabinets and Multirisk
- Containment and cans
- File cabinets
- Anti-fire equipments
- Showers and first-aid equipments