



TABLE / TROLLEY CABINET

P/N	Dimension H x W x D (mm)	Weight (kg)	Retention (L)
TA06	800 x 615 x 630	39	13
TA09	800 x 915 x 780	51	28
TA12	800 x 1215 x 780	78	38
TA15	800 x 1520 x 780	86	48

- Construction in 15/10th steel, epoxy paint white RAL 9010.
- 1 shelf adjustable in height.
- 1 removable Containment sump at the bottom.
- 2 hinged doors with automatic closing and key.
- Standardized warning labels.
- 4 swivel castors, 2 of which are lockable.
- Ventilation hole.



TROLLEY

P/N	Dimension H x W x D (mm)	Weight (kg)
TR06	800 x 615 x 630	34
TR09	800 x 917 x 780	50
TR12	800 x 1217 x 780	58
TR15	800 x 1517 x 780	66

- Construction in 15/10th steel.
- White epoxy paint RAL 9010.
- 1 shelf adjustable in height (the depth of the shelf is equal to $\frac{3}{4}$ of the depth of the table).
- 4 swivel castors, 2 of which are lockable.



© ECOSAFE reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 09/2021

Fume hoods, filtering cabinets - Ventilation

Lithium-ion cabinets

Flammable cabinets

Corrosive cabinets

Toxic cabinets and Multirisk

Containment and cans

File cabinets
Anti-fire equipments

Showers and first-aid equipments