

# RANGE 13



## VENTILATED SAFETY CABINETS FOR REAGENTS AND TOXIC PRODUCTS

### CONFORMITY

- In compliance with the health regulation articles (R.5132-66 & R.R132-68).
- Large warning labels in compliance with ISO 3864,92/58/CEE European directive and NF X08.003 norm.

### ADVANTAGES

- Integrated electric fan.
- Large storage capacity.

### PASSIVE SECURITY

- Built in steel 10/10<sup>th</sup>.
- White epoxy paint RAL 9010.
- «Toxic» and «Noxious» large warning labels (to specify with your order).
- 3 containment shelves adjustable in height with rack system.
- Removable retention tank at the bottom of cabinet (100 mm height).

### ACTIVE SECURITY

- Integrated electric fan : 220-240 volts - 0,1 A - 18 W; flow 90 m<sup>3</sup>/h; IPX4 protection ; sound level 37,5 Lp dB(A) (3m).
- Double handling doors with 3 points key-locking system.

### OPTIONS

- Additional shelf.
- Additional containment shelf in PVC.
- Sliding drawer.
- «Poison» chests (see below)

### «POISON» CHESTS

Can be used as wall mounted cabinets to store dangerous products or a poison chest in your safety cabinet ref MF4.

- Built in steel 10/10<sup>th</sup>.
- White epoxy paint RAL 9010.
- «Toxic» warning label in compliance with ISO 3864,92/58/CEE European directive and NF X08.003 norm.
- Shelves adjustable in height: 2 units for the model MF1 and 4 for the model MF2
- Handling door with key-locking system.
- Possibility to fix on the wall.
- Possibility to connect 2 models



^ MF4 model with MF1 «poison» chest option



MF2 ^



^ MF1

P/N	Description	Outside dim. H x W x D (mm)	Inside dim. H x W x D (mm)	Storage capacity (L)	Containment capacity (L)	Number of shelves	Maximum load by shelf <sup>(1)</sup>	W. (kg)
MF1	Poison chest 12 L	380 x 470 x 205	370 x 460 x 170	12	-	2	10	10
MF2	Poison chest 19 L	580 x 470 x 205	570 x 460 x 170	19	-	4	10	15
MF4	Cabinet with ventilation	1975 x 1000 x 500	1920 x 980 x 480	250	60	3	50	100
MF4SV	Cabinet without ventilation	1975 x 1000 x 500	1920 x 980 x 480	250	60	3	50	100
ESMF4	Additional shelf	-	-	-	-	-	-	6
SPMF4	Additional shelf in PVC	-	-	-	-	-	-	5
TCMF4	Sliding drawer	-	-	-	-	-	-	10

<sup>(1)</sup> Evenly distributed

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 08/2015

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk

Pesticide range

Fume hoods, filtering cabinets - Ventilation

Containment and Cans

Anti-fire equipments and Files cabinets

Showers and first-aid equipments