

## RANGE 8.B



# CORROSION RESISTANT CABINETS FOR ACIDS & BASES



▲ APC222

Secondary compartment



New slide rails with drawer stop



▲ APC62



▲ APC92



### CERTIFICATION

- Large warning labels in compliance with ISO 3864, 92/58/CEE Euro- pean directive and NF X08.003 norm.
- In compliance with recommendations of The National Institute for research and Safety at work (I.N.R.S.).

### ADVANTAGES

- Especially designed for acids and bases storage.
- Built in melamine-based wood with a great resistance to corrosive vapors aggressions.
- Cabinets with sliding drawers for acids or/and bases storage.

### PASSIVE SECURITY

- Built entirely in white melamine-based wood.
- «Corrosive», «Base» and «Toxic» large warning labels (only the 3 for APC107 model) .
- Adjustable retention sliding drawers for each major compartment.
- Independant compartments (except APC32 model) for acid, base or/ and toxic products storage (only for APC107 model).
- Each compartment is divided in a major and a secondary storage area (except APC62, APC112 and APC222 models).

### ACTIVE SECURITY

- High and low opening vent Ø 100 mm with connecting piece for possible exhaust.
- Door(s) with lock:
  - Handling door(s) for APC32, APC92 and APC107 models.
  - Sliding doors for APC62, APC112 and APC222 models.
- Adjustable leveling feet for enhance stability.

### Recycling AirBox Labopur - H50

Take advantage of **Labopur's** technology with any safety cabinet. Active coal filters absorbtion's capacity and filtration system are tested and certified NF X 15-211.



| P/N    | Outside Dim.<br>H x W x D (mm) | Number of compartment |           | Number of retention |                    | Storage capacity (kg) | Containment capacity (L) | Maximum load by shelf <sup>(1)</sup> (kg) | W. (kg) |
|--------|--------------------------------|-----------------------|-----------|---------------------|--------------------|-----------------------|--------------------------|---|---------|
|        |                                | Major                 | Secondary | Drawers in major    | Shelves in second. |                       |                          |   |         |
| APC32  | 670 x 635 x 615                | 1                     | 1         | 2                   | 2                  | 68                    | 28                       | 20  | 30      |
| APC62  | 670 x 1135 x 615               | 2                     | -         | 4                   | -                  | 96                    | 44                       | 20  | 60      |
| APC92  | 1850 x 635 x 615               | 2                     | 2         | 6                   | 6                  | 150                   | 84                       | 20  | 90      |
| APC112 | 900 x 1135 x 615               | 2                     | -         | 6                   | -                  | 144                   | 66                       | 20  | 110     |
| APC107 | 1850 x 635 x 615               | 3                     | 3         | 5                   | 5                  | 170                   | 70                       | 20  | 83      |
| APC222 | 1850 x 1135 x 615              | 2                     | -         | 8                   | -                  | 192                   | 88                       |   | 155     |
| E8A2   | Shelf / base trough sliding    |                       |           |                     |                    | 20                    | 11                       |   | 6       |

**N.B. :** The melamine-based wood may be affected by corrosive liquids attacks. Acids and bases bottles must be stored in the HDPE base trough. Moreover, all spills must be quickly cleaned. <sup>(1)</sup> Evenly distributed

The manufacturer reserves the right to discontinue or change specifications or designs at any time without notice and without incurring obligations. 05/2017

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk

Pesticide range

Fume hoods, filtering cabinets - Ventilation

Containment and Cans

Anti-fire equipments and Files cabinets

Showers and first-aid equipments