

1633R



Sicherheitsflasche - 1 L

ÜBEREINSTIMMUNG

ADVANTAGES

- Self-closing doors
- Retention equipment
- Storage of flammables and toxics

SONDERN AUCH...

- Fluchtabwesenheit, als das Fläschchen umstürzt
- Sealing cap with nitrile rubber joint
- Available in 2 colors (red or yellow)
- Available in different sizes and capacities
- Accurate dosing
- Double opening and sealing O-ring and spring
- Automatischer Verschluss/Schliessung
- Lid opening by finger pressure



TECHNISCHE DATEN

Referenz		1633R		Codenummer		73102910	
Referenz	Äußere Abmessungen HxBxT (mm)	Innere Abmessungen HxBxT (mm)	Lagerkapazität (L)	Auffangwanne	Regale	Tragfähigkeit	Gewicht (Kg)
1633R	280 x 80 x 80		1		0	0	0.20