



Regal aus Stahl - 4 Regale

ÜBEREINSTIMMUNG

ADVANTAGES

- Retention equipment
- Konform Artikel 2 Februar 1998
- Storage of flammables, corrosives and toxics

SONDERN AUCH...

- Zur Lagerung von korrosiven, giftigen und umweltschädlichen Flüssigkeiten
- Easy detection of leaks
- Containment sumps with edges
- Immediate view of the products
- High loading capacity
- Available in different sizes and storage capacity
- Allows the separation of different types of chemicals



TECHNISCHE DATEN

Referenz

RRA414

Codenummer

940310

Referenz	Äußere Abmessungen HxBxT (mm)	Innere Abmessungen HxBxT (mm)	Lagerkapazität (L)	Auffangwanne	Regale	Tragfähigkeit	Gewicht (Kg)
RRA414	2000 x 1030 x 390		144		4	0	45.00