

TRAC



NFX
15-211



Pack of 10 reactive tubes for corrosive vapors

ÜBEREINSTIMMUNG

ADVANTAGES

- Convenient for the measure of the filter's saturation
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics



TECHNISCHE DATEN

Referenz	TRAC	Codenummer	32041600
Referenz	Äußere Abmessungen HxBxT (mm)	Innere Abmessungen HxBxT (mm)	Gewicht (Kg)
TRAC			1.00



TRIONYX
DEPUIS 1971

138 avenue Henri BARBUSSE 92140 Clamart (France)
Tél. 01 46 45 80 00 - Fax : 01 46 45 85 98
www.trionyx.fr - trionyx[at]trionyx.fr