



Secured cabinet for the distribution of masks and hydroalcoholic gel - 20 lockers

COMPLIANCE

ADVANTAGES

- Storage of hydroalcoholic solution
- Storage of masks

BUT ALSO...

- Possibility to fix on the wall
- Individual key closing
- Individual label holder for each locker
- Individual reinforced door with noise damper



TECHNICAL FEATURES

Model
Cabinet materials
Doors color
Number of doors

ASC20C19
Robust steel and epoxy paint
Blue RAL 7035
20

Doors type
Hinges
Door closure
Number of compartments

Swinging
Right
Manual
20

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ASC20C19	1100 x 800 x 150	20x (210 x 193 x 121)	0		0	0	27.00