



## Secured cabinet for the distribution of masks and hydroalcoholic gel - 20 lockers

### COMPLIANCE

### ADVANTAGES

- Storage of hydroalcoholic solution
- Storage of masks

### BUT ALSO...

- Possibility to fix on the wall
- Individual key closing
- Individual label holder for each locker
- Individual reinforced door with noise damper



### TECHNICAL FEATURES

Model  
Cabinet materials  
Doors color  
Number of doors

**ASC20C19**  
**Robust steel and epoxy paint**  
**Blue RAL 7035**  
**20**

Doors type  
Hinges  
Door closure  
Number of compartments

**Swinging**  
**Right**  
**Manual**  
**20**

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ASC20C19	1100 x 800 x 150	20x (210 x 193 x 121)	0		0	0	27.00