



## Secured cabinet for the distribution of masks and hydroalcoholic gel – 35 lockers

### COMPLIANCE

### ADVANTAGES

- Storage of hydroalcoholic solution
- Storage of masks

### BUT ALSO...

- Possibility to fix on the wall
- Individual key closing
- Individual label holder for each locker
- Individual reinforced door with noise damper



### TECHNICAL FEATURES

Model	<b>ASC35C19</b>	Doors type	<b>Swinging</b>
Cabinet materials	<b>Robust steel and epoxy paint</b>	Hinges	<b>Right</b>
Doors color	<b>Blue RAL 7035</b>	Door closure	<b>Manual</b>
Number of doors	<b>35</b>	Number of compartments	<b>35</b>

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ASC35C19	1100 x 1380 x 150	35x (210 x 193 x 121)	0		0	0	43.00