



Secured cabinet for the distribution of masks and hydroalcoholic gel – 35 lockers

COMPLIANCE

ADVANTAGES

- Storage of hydroalcoholic solution
- Storage of masks

BUT ALSO...

- Possibility to fix on the wall
- Individual key closing
- Individual label holder for each locker
- Individual reinforced door with noise damper



TECHNICAL FEATURES

Model	ASC35C19	Doors type	Swinging
Cabinet materials	Robust steel and epoxy paint	Hinges	Right
Doors color	Blue RAL 7035	Door closure	Manual
Number of doors	35	Number of compartments	35

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ASC35C19	1100 x 1380 x 150	35x (210 x 193 x 121)	0		0	0	43.00