

AL157

# 1- door working cover style sheet steel safety cabinet

## COMPLIANCE

## ADVANTAGES

- Affordable equipment to store your phytosanitary products
- Compliant to the article 4 of the decree n°87-361 of 27 mai 1987
- Several levels of storage with retention
- Holes for the connection to the ventilation
- Allows the storage of your phytosanitary products safely
- Normalized labelling provided with the equipment
- Doors with key locking system



## TECHNICAL FEATURES

Model	<b>AL157</b>	Shipping dimensions HxWxD (mm)	<b>2000 x 1200 x 800</b>
Cabinet materials	<b>8/10th epoxy painted steel structure</b>	Ventilation	<b>Both inlet and outlet</b>
Doors color	<b>grey RAL 7035</b>	Openings diameter (mm)	<b>100</b>
Number of doors	<b>1</b>	Number of compartments	<b>1</b>
Doors type	<b>swing doors</b>	Containment sump	<b>1</b>
Hinges	<b>right</b>	HS Code	<b>94032080</b>
Door closure	<b>manual closing</b>		

## ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

## OPTIONS

- Perforated cover plate for 1-door cabinet (Réf : CAL155)
- Containment sump for 1-door cabinet (Réf : ESAL155)
- Containment sump in PVC for 1-door cabinet (Réf : SEAL155)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AL157	1800 x 500 x 420	1730 x 480 x 390	130	1	3	30	41.00