

2- door ventilated safety cabinet for toxics

COMPLIANCE

ADVANTAGES

- Compliant to the Public Health Code
- Complies with risk prevention and safety recommendations
- Storage of toxics and irritating products

BUT ALSO...

- Integrated electric fan
- High storage capacity
- Complies with risk prevention and safety recommendations
- Normalized labelling provided with the equipment
- 3-point locking system



TECHNICAL FEATURES

Model	MF4	Door closure	manual closing
Cabinet materials	10/10th epoxy painted steel structure	Shipping dimensions HxWxD (mm)	2175 x 1200 x 800
Doors color	white RAL 9010	Ventilation	Both inlet and outlet
Number of doors	2	Number of compartments	1
Doors type	swing doors	Containment sump	1
Hinges	right & left	HS Code	94032080

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Retention shelf for cabinets MF4 - MF4SV (Réf : ESMF4)
- Poison chest - 12L (Réf : MF1)
- Poison chest - 19L (Réf : MF2)
- Retention shelf in PVC for cabinets MF4 - MF4SV (Réf : SPMF4)
- Sliding drawer for cabinets MF4 - MF4SV (Réf : TCMF4)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
MF4	1975 x 1000 x 500	1920 x 980 x 480	250	1	3	50	100.00