



Yellow can 1 L

COMPLIANCE

ADVANTAGES

- Self-closing doors
- Retention equipment
- Storage of flammables and toxics

BUT ALSO...

- No leak when the flask overturns
- Sealing cap with nitrile rubber joint
- Available in 2 colors (red or yellow)
- Available in different sizes and capacities
- Accurate dosing
- Double opening and sealing O-ring and spring
- Automatic closing
- Lid opening by finger pressure



TECHNICAL FEATURES

| Model | 1633J | | | HS Code | 73102910 | | |
|-------|--------------------------------|--------------------------------|----------------------|------------------|----------|---------------|-------------|
| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Storage capacity (L) | Containment sump | Shelves | Load capacity | Weight (kg) |
| 1633J | 280 x 80 x 80 | | 1 | | 0 | 0 | 0.20 |