



Yellow can 1 L

COMPLIANCE

ADVANTAGES

- Self-closing doors
- Retention equipment
- Storage of flammables and toxics

BUT ALSO...

- No leak when the flask overturns
- Sealing cap with nitrile rubber joint
- Available in 2 colors (red or yellow)
- Available in different sizes and capacities
- Accurate dosing
- Double opening and sealing O-ring and spring
- Automatic closing
- Lid opening by finger pressure



TECHNICAL FEATURES

Model	1633J			HS Code	73102910		
Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
1633J	280 x 80 x 80		1		0	0	0.20