

AS5DD



2- door tall medicine cabinet

COMPLIANCE

ADVANTAGES

- Locking system
- Storage of pharmaceuticals products
- Storage of first-aid equipments



TECHNICAL FEATURES

Model	AS5DD	Door closure	manual closing
Cabinet materials	10/10th epoxy painted steel structure	Shipping dimensions HxWxD (mm)	2025 x 1200 x 800
Doors color	white RAL 9010	Ventilation	Inlet located at the bottom of the cabinet, outlet located at the top
Number of doors	2	Openings diameter (mm)	100
Doors type	swing doors	Number of compartments	1
Hinges	right & left	HS Code	94032080

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Eye-washer flow (L/min)	Shower flow (L/min)	Weight (kg)
AS5DD	1825 x 1100 x 545	1795 x 1095 x 480	0	0	155.00