

AS5DD



## 2- door tall medicine cabinet

### COMPLIANCE

### ADVANTAGES

- Locking system
- Storage of pharmaceuticals products
- Storage of first-aid equipments



### TECHNICAL FEATURES

Model	<b>AS5DD</b>	Door closure	<b>manual closing</b>
Cabinet materials	<b>10/10th epoxy painted steel structure</b>	Shipping dimensions HxWxD (mm)	<b>2025 x 1200 x 800</b>
Doors color	<b>white RAL 9010</b>	Ventilation	<b>Inlet located at the bottom of the cabinet, outlet located at the top</b>
Number of doors	<b>2</b>	Openings diameter (mm)	<b>100</b>
Doors type	<b>swing doors</b>	Number of compartments	<b>1</b>
Hinges	<b>right &amp; left</b>	HS Code	<b>94032080</b>

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Eye-washer flow (L/min)	Shower flow (L/min)	Weight (kg)
AS5DD	1825 x 1100 x 545	1795 x 1095 x 480	0	0	155.00