

MF15



## Wall-mounted medicine cabinet - 2 shelves

### COMPLIANCE

### ADVANTAGES

- Locking system
- Storage of pharmaceuticals products
- Storage of first-aid equipments



### TECHNICAL FEATURES

Model	<b>MF15</b>	Hinges	<b>right</b>
Cabinet materials	<b>10/10th epoxy painted steel structure</b>	Door closure	<b>manual closing</b>
Doors color	<b>white RAL 9010</b>	Shipping dimensions HxWxD (mm)	<b>580 x 600 x 800</b>
Number of doors	<b>1</b>	Number of compartments	<b>1</b>
Doors type	<b>swing doors</b>	HS Code	<b>94032080</b>

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Eye-washer flow (L/min)	Shower flow (L/min)	Weight (kg)
MF15	380 x 470 x 205	370 x 460 x 170	0	0	10.00