

MF15



Wall-mounted medicine cabinet - 2 shelves

COMPLIANCE

ADVANTAGES

- Locking system
- Storage of pharmaceuticals products
- Storage of first-aid equipments



TECHNICAL FEATURES

| | | | |
|-------------------|--|--------------------------------|------------------------|
| Model | MF15 | Hinges | right |
| Cabinet materials | 10/10th epoxy painted steel structure | Door closure | manual closing |
| Doors color | white RAL 9010 | Shipping dimensions HxWxD (mm) | 580 x 600 x 800 |
| Number of doors | 1 | Number of compartments | 1 |
| Doors type | swing doors | HS Code | 94032080 |

| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Eye-washer flow (L/min) | Shower flow (L/min) | Weight (kg) |
|-------|--------------------------------|--------------------------------|-------------------------|---------------------|-------------|
| MF15 | 380 x 470 x 205 | 370 x 460 x 170 | 0 | 0 | 10.00 |