

MF15



Wall-mounted medicine cabinet - 2 shelves

COMPLIANCE

ADVANTAGES

- Locking system
- Storage of pharmaceuticals products
- Storage of first-aid equipments



TECHNICAL FEATURES

Model	MF15	Hinges	right
Cabinet materials	10/10th epoxy painted steel structure	Door closure	manual closing
Doors color	white RAL 9010	Shipping dimensions HxWxD (mm)	580 x 600 x 800
Number of doors	1	Number of compartments	1
Doors type	swing doors	HS Code	94032080

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Eye-washer flow (L/min)	Shower flow (L/min)	Weight (kg)
MF15	380 x 470 x 205	370 x 460 x 170	0	0	10.00