

S1F

ISO
3864

EN
15154-1

EN
15154-2



Belt bag with 1 bottle 200 mL of sterile sodium chloride solution

COMPLIANCE

ADVANTAGES

- ISO 3864 compliant
- EN 15154-1 and EN 15154-2 compliant
- First-aid in case of chemical projection

BUT ALSO...

- Allows to rinse itself the eyes in case of chemical projection



TECHNICAL FEATURES

| Model | S1F | | HS Code | 28530010 | |
|-------|-----------------------------------|-----------------------------------|-------------------------|---------------------|-------------|
| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Eye-washer flow (L/min) | Shower flow (L/min) | Weight (kg) |
| S1F | | | 0 | 0 | 0.50 |