

S1F



Belt bag with 1 bottle 200 mL of sterile sodium chloride solution

COMPLIANCE

ADVANTAGES

- ISO 3864 compliant
- EN 15154-1 and EN 15154-2 compliant
- First-aid in case of chemical projection

BUT ALSO...

- Allows to rinse itself the eyes in case of chemical projection



TECHNICAL FEATURES

Model	S1F		HS Code	28530010	
Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Eye-washer flow (L/min)	Shower flow (L/min)	Weight (kg)
S1F			0	0	0.50