

Eye-wash kit delivered with 500 mL bottle with sodium chloride solution - 1 eye-wash bottle - 2 eye-pads - 1 scissors - 1 soft tissues (pack) - 1 adhe

COMPLIANCE

ADVANTAGES

- ISO 3864 compliant
- EN 15154-1 and EN 15154-2 compliant
- First-aid in case of chemical projection

BUT ALSO...

- Allows to rinse itself the eyes in case of chemical projection



TECHNICAL FEATURES

Model	KLY		HS Code	28530010	
Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Eye-washer flow (L/min)	Shower flow (L/min)	Weight (kg)
KLY			0	0	1.00