

Armadio di laboratorio anticorrosione - 30 L

CONFORMIDAD

Préconisations INRS

ADVANTAGES

- Made in France
- Conforme a las recomendaciones del INRS
- Storage of corrosives

SINO TAMBIÉN...

- PP containment sumps on sliding rails
- Conforme a las recomendaciones del INRS
- Built in cellular PVC with a great resistance to corrosive vapors and liquids
- Holes for the connection to the ventilation
- Especially designed to store acids and bases



CARACTERISTICAS TECNICAS

Modelo	APC31+	Cierre de las puertas	chiusura manuale
Conformidad	Préconisations INRS	Dimensiones de envío alt x anc x pro (mm)	870 x 600 x 800
Materiales del armario	PVC	Ventilación	orificio di ventilazione alto e bassi
Color de las puertas	bianco RAL 9010	Diámetro de ventilación (mm)	100
Número de puertas	1	Número de compartimentos	1
Tipo de puertas	porte battenti	Código Aduanero	94037000
Bisagras	cerniere a destra		

COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Cassettone di filtrazione acciaio Filtra CORG (Réf : CDF-ACORG)
- Cassettone di filtrazione PVC cellulare (Réf :)
- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Kit di raccordo per cassettone (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

OPTIONS

- Etagere/bac di ritenzione supplementare per gli scompartimenti acidi, basi, tossici (Réf : E8A2)

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
APC31+	630 x 635 x 578		40		0	20	30.00