

Armadio di laboratorio anticorrosione - 106 L

CONFORMIDAD

Préconisations INRS

ADVANTAGES

- Made in France
- Conforme a las recomendaciones del INRS
- Storage of corrosives

SINO TAMBIÉN...

- Conforme a las recomendaciones del INRS
- Especialmente diseñado to store acids and bases
- Built in melamine with a great resistance to corrosive vapors
- PP containment sumps on sliding rails
- Holes for the connection to the ventilation
- Cabinet with compartments to separate acids and bases



CARACTERISTICAS TECNICAS

Modelo
Conformidad

Materiales del armario
Color de las puertas
Número de puertas
Tipo de puertas
Bisagras

APC107+
Préconisations INRS

melamina
blanco RAL 9010
3
porte battenti
cerniere a sinistra

Cierre de las puertas
Dimensiones de envío alt x anc x pro (mm)
Ventilación
Diámetro de ventilación (mm)
Número de compartimentos
Código Aduanero

chiusura manuale
2100 x 1200 x 800
orifizio di ventilazione alto e bassi
100
3
94036000

COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Cassettone di filtrazione acciaio Filtra CORG (Réf : CDF-ACORG)
- Cassettone di filtrazione PVC cellulare (Réf :)
- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Gerla ed armadi hanno filtrazione - Ventilazione - Cassettone ha recirculation di aria senza filtro (Réf : H50C)
- Kit di raccordo per cassettone (Réf : KRC)

OPTIONS

- Etagere/bac di ritenzione supplementare per gli scompartimenti acidi, basi, tossici (Réf : E8A2)

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
APC107+	1900 x 635 x 615		108		0	20	83.00