

Armadio alta 2 portas

CONFORMIDAD

Articles R.5132-66 et R.5132-68 Code de la santé publique, Article 4 décret n°87-361

ADVANTAGES

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- Compliant to the article 4 of the decree n°87-361 of 27 mai 1987
- Several levels of storage with retention
- Holes for the connection to the ventilation
- Allows the storage of your phytosanitary products safely



CARACTERISTICAS TECNICAS

Modelo	AL307	Cierre de las puertas	chiusura manuale
Conformidad	Articles R.5132-66 et R.5132-68 Code de la santé publique, Article 4 décret n°87-361	Dimensiones de envío alt x anc x pro (mm)	2000 x 1200 x 800
Materiales del armario	acciaio 8/10	Ventilación	orifizio di ventilazione alto e bassi, laterale
Color de las puertas	grigio RAL 7035	Diámetro de ventilación (mm)	100
Número de puertas	2	Número de compartimentos	1
Tipo de puertas	porte battenti	Cubetas de retención	1
Bisagras	cerniere a destra e sinistra	Código Aduanero	94032080

COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Cassettone di filtrazione acciaio Filtra CORG (Réf : CDF-ACORG)
- Cassettone di filtrazione PVC cellulare (Réf :)
- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Gerla ed armadi hanno filtrazione - Ventilazione - Cassettone ha recirculation di aria senza filtro (Réf : H50C)
- Kit di raccordo per cassettoni (Réf : KRC)

OPTIONS

- Rejilla de Cubeta de retención - de armarios 2 puertas (Réf : CAL305)
- Cubeta de retención - de armarios 2 puertas (Réf : ESAL305)
- Cubeta de retención PVC - de armarios 2 puertas (Réf : SEAL305)

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
AL307	1800 x 920 x 420	1730 x 900 x 390	240	1	3	60	60.00