



# Armadio di sicurezza bicolore porta rosso - 155L

## CONFORMIDAD

Articles R.5132-66 et R.5132-68 Code de la santé publique

## ADVANTAGES

- Compliant to the Public Health Code
- Conforme a las recomendaciones del INRS
- Storage of flammables, corrosives and toxics

## SINO TAMBIÉN..

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- Conforme a las recomendaciones del INRS
- Allows the storage of your chemicals safely
- Affordable equipment to store your chemicals
- Top of the cabinet covered with a rubber carpet



## CARACTERISTICAS TECNICAS

Modelo	<b>AZ110R</b>	Cierre de las puertas	<b>chiusura manuale</b>
Conformidad	<b>Articles R.5132-66 et R.5132-68 Code de la santé publique</b>	Dimensiones de envío alt x anc x pro (mm)	<b>1200 x 1200 x 800</b>
Materiales del armario	<b>acciaio</b>	Ventilación	<b>orificio di ventilazione alto e bassi sulle porte</b>
Color de las puertas	<b>rosso RAL 3000</b>	Número de compartimentos	<b>1</b>
Número de puertas	<b>2</b>	Cubetas de retención	<b>1</b>
Tipo de puertas	<b>porte battenti</b>	Código Aduanero	<b>94032080</b>
Bisagras	<b>cerniere a destra e sinistra</b>		

## COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Opzione taglio copertura per raccordo (Réf : AZVENT)
- Cassettone di filtrazione acciaio Filtra CORG (Réf : CDF-ACORG)
- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Gerla ed armadi hanno filtrazione - Ventilazione - Cassettone ha recirculation di aria senza filtro (Réf : H50C)
- Kit di raccordo per cassettone (Réf : KRC)

## OPTIONS

- Carabottino per armadio AZ110 ed AZ300 (Réf : CE300)
- Mensola/bac per armadio AZ110 et AZ300 (Réf : ESE300)

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
AZ110R	1000 x 930 x 500	890 x 925 x 453	155	1	1	80	43.00