

MF4SV



Armadio ha reattivo e prodotti tossici senza ventilatore - 250 L

CONFORMIDAD

Articles R.5132-66 et R.5132-68 Code de la santé publique

ADVANTAGES

- Compliant to the Public Health Code
- Conforme a las recomendaciones del INRS
- Storage of toxics and irritating products

SINO TAMBIÉN...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- Conforme a las recomendaciones del INRS
- 3-point locking system



CARACTERISTICAS TECNICAS

Modelo	MF4SV	Cierre de las puertas	chiusura manuale
Conformidad	Articles R.5132-66 et R.5132-68 Code de la santé publique	Dimensiones de envío alt x anc x pro (mm)	2175 x 1200 x 800
Materiales del armario	acciaio 10/10	Ventilación	orificio di ventilazione sul tetto
Color de las puertas	bianco RAL 9010	Número de compartimentos	1
Número de puertas	2	Cubetas de retención	1
Tipo de puertas	porte battenti	Código Aduanero	94032080
Bisagras	cerniere a destra e sinistra		

COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Cassettone di filtrazione acciaio Filtra CORG (Réf : CDF-ACORG)
- Cassettone di filtrazione PVC cellulare (Réf : CDF-P)
- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Gerla ed armadi hanno filtrazione - Ventilazione - Cassettone ha recirculation di aria senza filtro (Réf : H50C)
- Kit di raccordo per cassettoni (Réf : KRC)

OPTIONS

- Mensola supplementare (Réf : ESMF4)
- Bagagliaio ha veleno (Réf : MF1)
- Bagagliaio ha veleno (Réf : MF2)
- Suretagere in polypropylene (Réf : SPMF4)
- Cassetto scorrevole (Réf : TCMF4)

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
MF4SV	1975 x 1000 x 500	1920 x 980 x 480	250	1	3	50	100.00