



# Armadio ventilato ha reattivo e prodotti tossici - 250 L

## CONFORMIDAD

Articles R.5132-66 et R.5132-68 Code de la santé publique

## ADVANTAGES

- Compliant to the Public Health Code
- Conforme a las recomendaciones del INRS
- Storage of toxics and irritating products

## SINO TAMBIÉN...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- Conforme a las recomendaciones del INRS
- 3-point locking system



## CARACTERISTICAS TECNICAS

Modelo	<b>MF4</b>	Cierre de las puertas	<b>chiusura manuale</b>
Conformidad	<b>Articles R.5132-66 et R.5132-68 Code de la santé publique</b>	Dimensiones de envío alt x anc x pro (mm)	<b>2175 x 1200 x 800</b>
Materiales del armario	<b>acciaio 10/10</b>	Ventilación	<b>Ventilatore elettrico integrato, orificio di ventilazione sull tetto</b>
Color de las puertas	<b>bianco RAL 9010</b>	Número de compartimentos	<b>1</b>
Número de puertas	<b>2</b>	Cubetas de retención	<b>1</b>
Tipo de puertas	<b>porte battenti</b>	Código Aduanero	<b>94032080</b>
Bisagras	<b>cerniere a destra e sinistra</b>		

## COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Kit di raccordo per cassettoni (Réf : KRC)

## OPTIONS

- Mensola supplementare (Réf : ESMF4)
- Bagagliaio ha veleno (Réf : MF1)
- Bagagliaio ha veleno (Réf : MF2)
- Suretagere in polypropylene (Réf : SPMF4)
- Cassetto scorrevole (Réf : TCMF4)

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
MF4	1975 x 1000 x 500	1920 x 980 x 480	250	1	3	50	100.00