

MR9+



## Armadio multirischi 90 min - EN 14470-1 - PVC cellulare

### CONFORMIDAD

EN 14470-1, EN 1363-1

### ADVANTAGES

- 105-minute fire resistant and made in France
- Storage of flammables, corrosives and toxics

### SINO TAMBIÉN...

- Built in cellular PVC for acids-bases-toxics compartments
- Holes for the connection to the ventilation
- Self-closing door if the temperature rises to 47°C
- Fixation for earthing (for flammable part)



### CARACTERISTICAS TECNICAS

|                        |   |   |  |
|------------------------|---|---|--|
| Modelo                 | <b>MR9+</b>                                   | Cierre de las puertas                     | <b>chiusura automatica e manuale</b>                 |
| Conformidad            | <b>EN 14470-1, EN 1363-1</b>                  | Dimensiones de envío alt x anc x pro (mm) | <b>2400 x 1300 x 800</b>                             |
| Materiales del armario | <b>melamina e acciaio 15/10 doppia parete</b> | Ventilación                               | <b>entrata dell'aria laterale e uscita sul tetto</b> |
| Color de las puertas   | <b>bianco RAL 9010</b>                        | Diámetro de ventilación (mm)              | <b>100</b>   |
| Número de puertas      | <b>4</b>                                      | Número de compartimentos                  | <b>4</b>   |
| Tipo de puertas        | <b>porte battenti</b>                         | Cubetas de retención                      | <b>1</b>   |
| Bisagras               | <b>cerniere a destra e sinistra</b>           | Código Aduanero                           | <b>94032080</b>                                      |

### COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Kit di raccordo per cassettoni (Réf : KRC)

### OPTIONS

- Etagere/bac di ritenzione supplementare per gli scompartimenti acidi, basi, tossici (Réf : E8A2)
- Carabottino per traghettino di ritenzione per armadio 1 porta - L 635 - P 615 mm (Réf : B148)
- Carabottino per traghettino di ritenzione per armadio 1 porta - L 635 - P 615 mm (Réf : C148)
- Mensola di ritenzione supplementare per lo scompartimento infiammabile (Réf : E48)
- Cassetto/traghettino scorrevole per armadio 1 porta - L 635 - P 615 mm (Réf : T148)
- Mensola di ritenzione supplementare per lo scompartimento infiammabile (Réf : E48)

| Modelo | Dimensiones exteriores alt x anc x pro (mm) | Dimensiones interiores alt x anc x pro (mm) | Capacidad de almacenamiento (L) | Cubetas de retención | Estanterías | Carga admisible | Peso (kg) |
|--------|---|---|---------------------------------|----------------------|-------------|-----------------|-----------|
| MR9+   | 2205 x 1270 x 615                           |   | 210                             | 1                    | 3           | 30              | 400.00    |