

MF25



## Armadio ha farmacia murale, grande modello,

### CONFORMIDAD

### ADVANTAGES

- Storage of pharmaceuticals products
- Storage of first-aid equipments



### CARACTERISTICAS TECNICAS

Modelo	<b>MF25</b>	Bisagras	<b>cerniere a destra</b>
Materiales del armario	<b>acciaio 10/10</b>	Cierre de las puertas	<b>chiusura manuale</b>
Color de las puertas	<b>bianco RAL 9010</b>	Dimensiones de envío alt x anc x pro (mm)	<b>780 x 600 x 800</b>
Número de puertas	<b>1</b>	Número de compartimentos	<b>1</b>
Tipo de puertas	<b>porte battenti</b>	Código Aduanero	<b>94032080</b>

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavajojos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
MF25	580 x 470 x 205	570 x 460 x 170	0	0	15.00