

MF15



Armadio ha farmacia murale, piccolo modello,

CONFORMIDAD

ADVANTAGES

- Storage of pharmaceuticals products
- Storage of first-aid equipments



CARACTERISTICAS TECNICAS

Modelo	MF15	Bisagras	cerniere a destra
Materiales del armario	acciaio 10/10	Cierre de las puertas	chiusura manuale
Color de las puertas	bianco RAL 9010	Dimensiones de envío alt x anc x pro (mm)	580 x 600 x 800
Número de puertas	1	Número de compartimentos	1
Tipo de puertas	porte battenti	Código Aduanero	94032080

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavajjos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
MF15	380 x 470 x 205	370 x 460 x 170	0	0	10.00