

MF15



## Armadio ha farmacia murale, piccolo modello,

### CONFORMIDAD

### ADVANTAGES

- Storage of pharmaceuticals products
- Storage of first-aid equipments



### CARACTERISTICAS TECNICAS

|                        |                        |                                           |                          |
|------------------------|------------------------|-------------------------------------------|--------------------------|
| Modelo                 | <b>MF15</b>            | Bisagras                                  | <b>cerniere a destra</b> |
| Materiales del armario | <b>acciaio 10/10</b>   | Cierre de las puertas                     | <b>chiusura manuale</b>  |
| Color de las puertas   | <b>bianco RAL 9010</b> | Dimensiones de envío alt x anc x pro (mm) | <b>580 x 600 x 800</b>   |
| Número de puertas      | <b>1</b>               | Número de compartimentos                  | <b>1</b>                 |
| Tipo de puertas        | <b>porte battenti</b>  | Código Aduanero                           | <b>94032080</b>          |

| Modelo | Dimensiones exteriores alt x anc x pro (mm) | Dimensiones interiores alt x anc x pro (mm) | Caudal de lavajjos (L/min) | Caudal de la ducha de seguridad (L/min) | Peso (kg) |
|--------|---------------------------------------------|---------------------------------------------|----------------------------|-----------------------------------------|-----------|
| MF15   | 380 x 470 x 205                             | 370 x 460 x 170                             | 0                          | 0                                       | 10.00     |