

SM1F

ISO
3864EN
15154-1EN
15154-2

Stazione murale 1 flacone di 500 ml per gli occhi

CONFORMIDAD

ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection

SINO TAMBIÉN...

- Allows to rinse itself the eyes in case of chemical projection



CARACTERISTICAS TECNICAS

Modelo

SM1F

Código Aduanero

28530010

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavajojos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
SM1F			0	0	0.90

**TRIONYX**
DEPUIS 1971

138 avenue Henri BARBUSSE 92140 Clamart (France)
Tél. 01 46 45 80 00 - Fax : 01 46 45 85 98
www.trionyx.fr - trionyx[at]trionyx.fr