

SM1F

ISO
3864EN
15154-1EN
15154-2

Stazione murale 1 flacone di 500 ml per gli occhi

CONFORMIDAD

ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection

SINO TAMBIÉN...

- Allows to rinse itself the eyes in case of chemical projection



CARACTERISTICAS TECNICAS

Modelo

SM1F

Código Aduanero

28530010

| Modelo | Dimensiones exteriores alt x anc x pro (mm) | Dimensiones interiores alt x anc x pro (mm) | Caudal de lavajojos (L/min) | Caudal de la ducha de seguridad (L/min) | Peso (kg) |
|--------|---|---|-----------------------------|---|-----------|
| SM1F | | | 0 | 0 | 0.90 |