

SM2F

ISO  
3864

EN  
15154-1

EN  
15154-2



## Stazione murali 2 flaconi di 500 ml per gli occhi

### CONFORMIDAD

### ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection

### SINO TAMBIÉN...

- Allows to rinse itself the eyes in case of chemical projection



### CARACTERISTICAS TECNICAS

Modelo

SM2F

Código Aduanero

28530010

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavaojos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
SM2F			0	0	1.70