

SP2F



Stazione murali 2 flaconi: 500 ml per gli occhi e 200 ml per il corpo

CONFORMIDAD

ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection

SINO TAMBIÉN...

- Allows to rinse itself the eyes in case of chemical projection



CARACTERISTICAS TECNICAS

Modelo	SP2F		Código Aduanero	28530010	
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavaojos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
SP2F			0	0	1.20